



SmartCare Remote Management

Feature compatibility matrix

							
Feature	CONNEX 360	CONNEX Vital Signs Monitor	CONNEX Integrated Wall System	CONNEX SPOT Monitor	RETINAVUE 700 Imager	CENTRELLA Smart+ Bed and PROGRESSA Smart+ Bed	DYNAMO Series Stretchers
Actionable Insights	X					X	X
Remote log retrieval	X	X	X	X	X		
Preventive maintenance tracking	X	X	X	X		X	X
View error code notifications						X	X
Component health alerts	X	X	X	X			
Access and print reports	X	X	X	X	X	X	X
Access and deploy configuration files	X	X	X	X			
Access and deploy custom data files	X	X	X	X			
Upgrade Firmware	X			X	X	X	X
Locate connected Baxter products via access point (AP) mapping information	X	X	X	X		X	X

For information about any Baxter product, contact Baxter Technical Support: baxter.com/contact-us.

Baxter, Centrella, Connex, Dynamo Series, Progressa, and RetinaVue are trademarks of Baxter International, Inc. or its subsidiaries.
80030603 D, Revision date: 2026-02