


Scale-Tronix

Scale dimensions for Equipment Planners

2001 SlingScale [\(Catalog, page 16\)](#)



Weight 85

2002 SlingScale [\(Catalog ,page 14\)](#)



Weight
100lbs

2250CE Pod Scale


[\(Catalog, page 17\)](#)



Weight 25lbs

4302 Organ / Specimen Scale

(Catalog, page 28)



Weight 9Lbs.

4502 Diaper Scale

(catalog, page 25)



Weight 9Lbs.

4802D

Pediatric Scale / Infant Scale [\(catalog, page23\)](#)



48800-1 Cart for Infant Scale

Weight (Scale 19Lbs),(Cart 41Lbs).


5002 Stand-On Scale [\(Catalog, page 3\)](#)



Weight 62Lbs.


5102 Physician's Stand-On Scale

(Catalog,6-7)



5107 Bariatric Physician's Stand-On Scale

[\(catalog, page 9\)](#)



5122 Low Profile Stand-On Scale


[\(catalog,page7\)](#)



Weight 20Lbs.

5125 Portable Stand-on Scale


(catalog, page 6)



Weight 12Lbs.

5127 Bariatric Portable Stand-On Scale

(Catalog, page 8)



Weight 18lbs

5152 In-Floor Stand-On Scale

(Catalog page 21)




18"



Flush readout

11"


15"



Surface mounted readout

10.25"

14.125"



22.5"

21"

2.5"

Weight 50Lbs.

5202 Stow-A-Weigh Stand-on Scale

(catalog, page 4)



For Bracket Location
see next slide.



18"
Weight 37Lbs.


5202 Bracket Location

WALL BRACKET MOUNTING INSTRUCTIONS

14 INCH BRACKET P/N 600019

5202 STAND-ON

LOCATION OF BRACKET MOUNTING
HOLES WITH RESPECT TO FLOOR



5602 Mini Platform Stand-On Scale

[\(catalog, page 7\)](#)



Weight 14Lbs.

5702 Bariatric Stand-On Scale

(catalog, page 8)



Weight 85Lbs.

6102 Flush-Mounted In-Floor Scale (computer enhanced)

(Catalog, page 20)



36"

37"

32"

34"

Flush readout

11"



15"


or Surface mounted readout

10.25"



14.125"

3.5"



38"


Weight 115Lbs.

6102D Dual platform scale with optional handrail.

[\(catalog, page 21\)](#)




6154 Flush-Mounted In-Floor Scale (computer enhanced) ([catalog, page #](#))



Weight 166Lbs.

6002 Wheelchair Scale


(catalog, page 10)



Weight 80Lbs.

6202 Stow-A-Weigh Wheelchair scale

[\(catalog, page 12\)](#)




Weight 89Lbs.

6202 Bracket Location

14 INCH BRACKET P/N 600019

6202 WHEELCHAIR SCALES

LOCATION OF BRACKET MOUNTING
HOLES WITH RESPECT TO FLOOR



6202D Dual-Stow-A-Weigh Scale




See next slide for bracket locations

6202D Bracket location

DETAILS FOR 6202D DUAL SCALE WEIGHING SYSTEM USING 14 INCH BRACKETS

LOCATION OF BRACKET MOUNTING
HOLES WITH RESPECT TO FLOOR



6502 Chair Scale

(catalog, page 26)



Weight 94Lbs.

6702W Oversized Wheelchair Scales

(catalog, page 10)



Weight 115Lbs.

6702SP Oversized Wheelchair Scales

(catalog, page 10)



Weight 122Lbs.

6702 Oversized Wheelchair Scales

(catalog, page 10)



Weight 97Lbs.

6772 Bariatric Wheelchair Scale


([catalog, page 11](#))



Weight 55.5Lbs.

845010W Wall Mounted Height Gauge

([catalog, page 13](#))



Weight 4Lbs.