

Welch Allyn ProBP 2400

Serial Number Location & SAP Lookup

Important:

When searching for ProBP 2400 equipment records in SAP:

Enter asterisk * then enter last 8 digits/characters of SN as stated by customer. This will surface the associated ProBP2400 equipment record in SAP.

Special Note: There are 3,000 devices provided by manufacturer that when scanned at shipping put "JC121" in front of the true device serial number. We should not share w/ the customer their serial number begins with JC121.



*GTIN and SN Label
(21) Is NOT a part of the SN*