

CSM IWS, CSM 777, CSM 77710, 411445

Part # 411445 is used to mount the Connex Spot Monitor to the wallboard for the Integrated Diagnostic System.

If a customer wants to mount their CSM to their wall system board, sell them 411445

Per Product Manager, we may sell this direct (even tho it is not coded as DPA)

What is included when PN 411445 ships: CSM Mounting bracket and hardware kit to mount bracket to wall board, Wall Bin w/ SpO2 Spool and mounting hardware

*** Note the power supply ships with the CSM monitor that is sold separate.

Material Document General					
It...	ICt	Component	Component description	Quan...	Un
0010 L	716651		WF BIN MOUNTING BRACKET	1	EA
0020 L	411427		LEFT BIN, MACHINED	1	EA
0030 L	723812		VESA MOUNT, TILT, LCD-203T	1	EA
0040 L	723791		SPACER, 35W BIN	2	EA
0050 L	406941		SPO2 SPOOL KIT, ACM	1	EA
0080 L	719950		LABEL, BLANK 1" x 4", ROLL	1	EA
0090 L	411433		HARDWARE KIT, CSM WALL	1	EA
0100 L	703831		BUBBLEWRAP 1/2X24	1	EA
0110 L	712105		MODPG PUMP FOAM	2	EA
0120 L	M11039		3X4 BAG	1	EA
0130 L	M09504		OVER-PACK BOX	1	EA

Link to the **mounting instruction for images of the parts**- this is on website and may be shared with customers: [See document 80020744](#)

PART 411445 - VESA MOUNT FOR CSM FOR INTEGRATED WALL SYSTEM

Create Standard Order: Overview

Standard Order: Net value: 215.00 USD

Sold-To Party: 10101235 Template Med End User Credit Card / Skaneateles Falls NY 13153 Find

Ship-To Party: 10101235 Template Med End User Credit Card / Skaneateles Falls NY 13153

PO Number: PO date:

Sales Item overview Item detail Ordering party Procurement Shipping Reason for rejection

Req. delv.date: 07/30/2018 Deliver Plant:

☐ Complete div. Total Weight: 0 KG

Delivery block: Volume: 0.000

Billing block: Pricing date: 07/27/2018

Payment card: Exp.date: CCARD

Card Verif.Code: Entry mode:

Payment terms: 0098 Credit Card Sale Incoterms: CPT Factory of Origin

Order reason:

Main items

Item	Material	Order Quantity	Description	Un	ItCa	Prit	S	DGIP	HL	Item	D First date
10	411445		1 KIT_BIN, SPCR, BRKT, SPL, PS, VESA, HD, EA	EA	TAN	1000					D 08/03/2018
											D 07/30/2018
											D 07/30/2018
											D 07/30/2018
											D 07/30/2018